

# Bio Image Operation script operations (v1.6.440.1 / 2020-11-14)

**SetPath** (Path)

Set path for relative file paths (by default path of current script file)

- Path: File path ("path")

**CreateImage** (Width, Height, ColorMode, Red, Green, Blue)

Create a new image

- Width: Width in pixels (numeric value)
- Height: Height in pixels (numeric value)
- ColorMode: Color mode (GrayScale, Color, ColorAlpha)
- Red: Red color component (numeric value between 0 and 1)
- Green: Green color component (numeric value between 0 and 1)
- Blue: Blue color component (numeric value between 0 and 1)

**OpenImage** (Path, Start, Length, Interval)

Open image file(s) for processing, accepts file name pattern

- Path: File path ("path")
- Start: Start (time reference as (hours:)minutes:seconds, or frame number)
- Length: Length (time reference as (hours:)minutes:seconds, or frame number)
- Interval: Interval in number of frames (numeric value)

**OpenVideo** (Path, API, Start, Length, Interval)

Open video file(s) and process frames, accepts file name pattern (ffmpeg formats supported)

- Path: File path ("path")
- API: OpenCV API code (See OpenCV API codes) (numeric value)
- Start: Start (time reference as (hours:)minutes:seconds, or frame number)
- Length: Length (time reference as (hours:)minutes:seconds, or frame number)
- Interval: Interval in number of frames (numeric value)

**OpenCapture** (API, Path, Source, Width, Height, Interval)

Open capturing from video (IP) path or camera source

- API: OpenCV API code (See OpenCV API codes) (numeric value)
- Path: File path ("path")
- Source: Camera source (#) (numeric value)
- Width: Width in pixels (numeric value)
- Height: Height in pixels (numeric value)
- Interval: Interval in number of frames (numeric value)

**SaveImage** (Path, Label, Start, Length)

Save image to file

- Path: File path ("path")
- Label: Label id (string)

- **Start:** Start (time reference as (hours:)minutes:seconds, or frame number)
- **Length:** Length (time reference as (hours:)minutes:seconds, or frame number)

**SaveVideo** (Path, Label, Start, Length, Fps, Codec)

Create video file and save image to video file (supports installed encoders)

- **Path:** File path ("path")
- **Label:** Label id (string)
- **Start:** Start (time reference as (hours:)minutes:seconds, or frame number)
- **Length:** Length (time reference as (hours:)minutes:seconds, or frame number)
- **Fps:** Frames per second (numeric value)
- **Codec:** Video encoding codec (4 character codec reference (FOURCC))

**ShowImage** (Label, Display)

Show image on screen (low priority screen updates)

- **Label:** Label id (string)
- **Display:** Display id (number 1 - 4)

**StoreImage** (Label)

Store current image in memory

- **Label:** Label id (string)

**GetImage** (Label)

Get specified stored image from memory

- **Label:** Label id (string)

**Grayscale** (Label)

Convert image to gray scale

- **Label:** Label id (string)

**Color** (Label)

Convert image to color

- **Label:** Label id (string)

**ColorAlpha** (Label)

Convert image to color with alpha channel

- **Label:** Label id (string)

**GetSaturation** (Label)

### Extract saturation from image

- Label: Label id (string)

### GetHsValue (Label)

### Extract (HSV) Value from image

- Label: Label id (string)

### GetHsLightness (Label)

### Extract (HSL) Lightness from image

- Label: Label id (string)

### Scale (Width, Height, Label)

### Scale image

- Width: Width in pixels (numeric value)
- Height: Height in pixels (numeric value)
- Label: Label id (string)

### Crop (X, Y, Width, Height, Label)

### Crop image

- X: X position in pixels (numeric value)
- Y: Y position in pixels (numeric value)
- Width: Width in pixels (numeric value)
- Height: Height in pixels (numeric value)
- Label: Label id (string)

### Mask (Label)

### Perform mask on current image

- Label: Label id (string)

### Threshold (Label, Level)

### Convert image to binary using threshold level, or in case not provided using automatic Otsu method

- Label: Label id (string)
- Level: Threshold value (numeric value between 0 and 1)

### Difference (Label)

### Perform difference of current image and specified image

- Label: Label id (string)

### **DifferenceAbs (Label)**

Perform absolute difference of current image and specified image

- Label: Label id (string)

### **Add (Label)**

Adds specified image to current image

- Label: Label id (string)

### **Multiply (Factor)**

Perform multiplication of all color channels by specified factor

- Factor: Multiplication factor (numeric value)

### **Invert (Label)**

Invert image

- Label: Label id (string)

### **UpdateBackground (Label, Weight)**

Add image to the adaptive background buffer

- Label: Label id (string)
- Weight: Weight value (numeric value between 0 and 1)

### **UpdateAverage (Label, Weight)**

Add image to the average buffer

- Label: Label id (string)
- Weight: Weight value (numeric value between 0 and 1)

### **ClearSeries ()**

Clear image series buffer

### **AddSeries (Label, Maximum)**

Add image to image series buffer

- Label: Label id (string)
- Maximum: Maximum number of images to keep (numeric value)

### **GetSeriesMedian ()**

Retrieve median image of image series buffer

**AddAccum** (Label, AccumMode)

Add image to the accumulative buffer

- Label: Label id (string)
- AccumMode: Accumulation mode (Age, Usage)

**GetAccum** (Power, Palette)

Retrieve the accumulative buffer and convert to image

- Power: Exponential power of value range (1E-[power] ... 1) (numeric value)
- Palette: Palette (GrayScale, Heat, Rainbow)

**CreateClusters** (Tracker, MinArea, MaxArea)

Create clusters; auto calibrate using initial images if no parameters specified

- Tracker: Tracker id (string)
- MinArea: Minimum area in number of pixels (numeric value)
- MaxArea: Maximum area in number of pixels (numeric value)

**CreateTracks** (Tracker, MaxMove, MinActive, MaxInactive)

Create cluster tracking; auto calibrate using initial images if no parameters specified

- Tracker: Tracker id (string)
- MaxMove: Maximum movement distance (single frame) (numeric value)
- MinActive: Minimum number of frames being active before state is active (numeric value)
- MaxInactive: Maximum number of frames being inactive before state is inactive (numeric value)

**CreatePaths** (Tracker, Distance)

Create common path usage

- Tracker: Tracker id (string)
- Distance: Maximum path distance (numeric value)

**DrawClusters** (Label, Tracker, DrawMode)

Draw clusters

- Label: Label id (string)
- Tracker: Tracker id (string)
- DrawMode: (Combination of) draw mode(s) (None, Point, Circle, Box, Angle, Label, Labeln, Track, Tracks, Fill, ClusterDefault, TracksDefault)

**DrawTracks** (Label, Tracker, DrawMode)

Draw tracked clusters

- Label: Label id (string)

- Tracker: Tracker id (string)
- DrawMode: (Combination of) draw mode(s) (None, Point, Circle, Box, Angle, Label, Labeln, Track, Tracks, Fill, ClusterDefault, TracksDefault)

**DrawPaths** (Label, Tracker, PathDrawMode, Power, Palette)

Draw common paths

- Label: Label id (string)
- Tracker: Tracker id (string)
- PathDrawMode: Path draw mode (Age, Usage, Usage2, Links, LinksMove)
- Power: Exponential power of value range (1E-[power] ... 1) (numeric value)
- Palette: Palette (GrayScale, Heat, Rainbow)

**SaveClusters** (Path, Tracker, Format, Contour)

Save clusters to CSV file

- Path: File path ("path")
- Tracker: Tracker id (string)
- Format: Output format (ByTime, ByLabel, Split)
- Contour: Extract contours (true / false)

**SaveTracks** (Path, Tracker, Format, Contour)

Save cluster tracking to CSV file

- Path: File path ("path")
- Tracker: Tracker id (string)
- Format: Output format (ByTime, ByLabel, Split)
- Contour: Extract contours (true / false)

**SavePaths** (Path, Tracker)

Save paths to CSV file

- Path: File path ("path")
- Tracker: Tracker id (string)

**ShowTrackInfo** (Tracker, Display)

Show tracking information on screen

- Tracker: Tracker id (string)
- Display: Display id (number 1 - 4)

**DrawTrackInfo** (Label, Tracker)

Draw tracking stats on image

- Label: Label id (string)
- Tracker: Tracker id (string)

**SaveTrackInfo** (Path, Tracker)

### Save tracking information to CSV file

- Path: File path ("path")
- Tracker: Tracker id (string)

### SaveTrackLog (Path, Tracker)

### Save tracking log to CSV file

- Path: File path ("path")
- Tracker: Tracker id (string)

### DrawLegend (Label, Display, Position)

### Draw legend

- Label: Label id (string)
- Display: Display id (number 1 - 4)
- Position: Draw position (Full, TopLeft, BottomLeft, TopRight, BottomRight)

### Wait (MS)

### Pause execution for a period (1000 ms default)

- MS: Time in milliseconds (numeric value)

### Debug ()

### Debug mode

**Arguments:**    Required        Optional