

# BioImageOperation script operations

**SetPath**(Path)

- Set path for relative file paths (by default path of current script file)

**CreateImage**(Width, Height, ColorMode, Red, Green, Blue)

- Create a new image

**OpenImage**(Path, Start, Length, Interval)

- Open image file(s) for processing, accepts file name pattern, start/length #

**OpenVideo**(Path, API, Start, Length, Interval)

- Open video file(s) and process frames, accepts file name pattern, start/length frames or time (ffmpeg formats supported) (API: See OpenCV API codes)

**OpenCapture**(API, Path, Source, Width, Height, Interval)

- Open capturing from video (IP) path or camera source (#) (API: See OpenCV API codes)

**SaveImage**(Path, Label, Start, Length)

- Save image to file, start/length frames or time

**SaveVideo**(Path, Label, Start, Length, Fps, Codec)

- Create video file and save image to video file, start/length frames or time, using fourcc codec string (supports installed encoders)

**ShowImage**(Label, Display)

- Show image on screen (low priority screen updates) (display: 1 - 4)

**StoreImage**(Label)

- Store current image in memory

**GetImage**(Label)

- Get specified stored image from memory

**Grayscale**(Label)

- Convert image to gray scale

**Color**(Label)

- Convert image to color

**ColorAlpha**(Label)

- Convert image to color with alpha channel

**GetSaturation**(Label)

- Extract saturation from image

**GetHsValue**(Label)

- Extract (HSV) Value from image

**GetHsLightness**(Label)

- Extract (HSL) Lightness from image

**Scale**(Width, Height, Label)

- Scale image

**Crop**(X, Y, Width, Height, Label)

- Crop image

**Mask**(Label)

- Perform mask on current image

**Threshold**(Label, Level)

- Convert image to binary using threshold level, or in case not provided using automatic Otsu method

**Difference**(Label)

- Perform difference of current image and specified image

**DifferenceAbs**(Label)

- Perform absolute difference of current image and specified image

**Add**(Label)

- Adds specified image to current image

**Multiply**(Factor)

- Perform multiplication of all color channels by specified factor

**Invert**(Label)  
•Invert image

**UpdateBackground**(Label, Weight)  
•Add image to the adaptive background buffer

**UpdateAverage**(Label, Weight)  
•Add image to the average buffer

**AddSeries**(Label, Maximum)  
•Add image to image series buffer

**GetSeriesMedian**()  
•Retrieve median image of image series buffer

**AddAccum**(Label, AccumMode)  
•Add image to the accumulative buffer (AccumMode: Age, Usage)

**GetAccum**(Power, Palette)  
•Retrieve the accumulative buffer and convert to image (Palette: Grayscale, Heat, Rainbow)

**CreateClusters**(Tracker, MinArea, MaxArea)  
•Create clusters; auto calibrate using initial images if no parameters specified

**CreateTracks**(Tracker, MaxMove, MinActive, MaxInactive)  
•Create cluster tracking; auto calibrate using initial images if no parameters specified

**CreatePaths**(Tracker, Distance)  
•Create common path usage

**DrawClusters**(Label, Tracker, DrawMode)  
•Draw clusters (DrawMode: Point|Circle|Box|Angle|Label|LabelIn|Fill)

**DrawTracks**(Label, Tracker, DrawMode)  
•Draw tracked clusters (DrawMode: Point|Circle|Box|Angle|Label|LabelIn|Track|Tracks)

**DrawPaths**(Label, Tracker, PathDrawMode, Power, Palette)  
•Draw common paths (PathDrawMode: Age, Usage, Links, LinksMove) (Palette: Grayscale, Heat, Rainbow)

**SaveClusters**(Path, Tracker, ByLabel)  
•Save clusters to file (CSV format)

**SaveTracks**(Path, Tracker, ByLabel)  
•Save cluster tracking to file (CSV format)

**SavePaths**(Path, Tracker)  
•Save paths to file (CSV format)

**ShowTrackInfo**(Tracker, Display)  
•Show tracking information on screen (Display: 1 - 4)

**DrawTrackInfo**(Label, Tracker)  
•Draw tracking stats on image

**SaveTrackInfo**(Path, Tracker)  
•Save tracking information to file (CSV format)

**SaveTrackLog**(Path, Tracker)  
•Save tracking log to file (CSV format)

**DrawLegend**(Label, Display, Position)  
•Draw legend (Display: 1 - 4 or Position on image: TopLeft, BottomLeft, TopRight, BottomRight)

**Wait**(MS)  
•Pause execution for a period (1000 ms default)

**Debug**()  
•Debug mode

**Arguments:**      Required      Optional